
INdEX columns

Guide to Safe Operation

Contents

Dansk	5
Deutsch	7
English	9
Español	11
Français	13
Italiano	15
Nederlands	17
Portuguès	19
Suomalainen	21
Svensk	23
Decontamination Report	25

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
INdEX COLUMNS Guide to Safe Operation

1 Introduction

This guide contains essential information which should be observed during installation, operation and maintenance. It is imperative that the User Manual is read by the operator before assembly, operation and maintenance. The Guide to Safe Operation and the User Manual should always be on hand whenever the column is used.

2 Labels relating to safety precautions

There is a risk of personal injury or damage to the column if the safety precautions in this document are not followed. Warning labels supplied with the column are shown to the right. An extra package of labels supplied with the column contains the warning text in Danish, Dutch, Finnish, French, German, Italian, Portuguese, Spanish and Swedish. Before using the column, ensure that the label with the relevant language for the country concerned is affixed to the column.

3 Safety precautions for maintenance, inspection and installation work

The end user must ensure that all maintenance, inspection and installation work is carried out by qualified personnel who are adequately trained and understand the operating instructions.

4 Specific points for the use of GE Healthcare's INdEX columns

INdEX™ columns should only be used for the purpose for which they are sold. The design pressure of 3 bar should never be exceeded. The working temperature is from +4 to +40 °C.

4.1 Chemical resistance

The table is intended as a guide to the resistance of the materials in INdEX columns to chemicals used as solvents. The information has been compiled from published material from several sources, not from individual tests on the components of the column. The operator should therefore use this list only as a guide to the level of chemical resistance exhibited by each material.

In general, the use of the following chemicals should be avoided:

- Extreme oxidizers (such as peroxides)
- Fluorine and halogenated compounds
- Chlorinated solvents (such as methylene chloride)
- Esters (such as acetates)
- Aromatic hydrocarbons (such as toluene)
- High concentration of strong acids

4.2 Emptying the column

Always make sure that there is enough free space for the adaptor to rise to its full extent when emptying the column, otherwise there is a risk that the adaptor will jam/stick or that the design pressure will be exceeded. To avoid exceeding the design pressure, install a pressure relief valve. For positional installation of the pressure relief valve refer to Figure 3, Schematic arrangement of the column, valves and pump, in the User Manual.

Chemical resistance of materials in INdEX columns.

Substance	Concentration by volume	60–90 days ¹
Acetic acid	10%	See note ²
Acetonitrile	5%	OK
Acetone	10%	OK
Ethanol	70%	OK
Ethylene glycol	50%	OK
Glycerol	100%	OK
Hydrochloric acid	0.1 M	See note ^{2,3}
Methanol	70%	OK
Nitric acid	0.1 M	See note ²
n-Propanol	70%	OK
Sodium chloride	2 M	See note ⁴
Sodium hydroxide	2 M	See note ⁵
Trifluoroacetic acid	0.1%	OK
Triton X-100	100%	OK
UREA	8 M	OK

¹ The test does not include PVC tubing.

² Polyamide nets cannot be used.

³ Not longer than 4 hours.

⁴ Can be used under normal running conditions. Do not use NaCl in storage solutions. Please note that NaCl can cause corrosion on stainless steel in acid solutions, pH<4.

⁵ Not longer than 3 hours.

4.3 Assembly recommendations

- Do not place fingers under the lid when mounting it onto the column as there is a risk of personal injury.
- If the tubing on the hydraulic inlet and outlet valves is incorrectly assembled, pressurised liquid may spurt out. To avoid this, refer to the figures in the User Manual, Chapter 5.
- Ensure that the tubing is long enough to allow for complete movement of the adaptor. If tubing is too short, it can be torn off the column, releasing liquid that may not be environmentally friendly.
- After reassemble the column ensure that the column is tight and not leaking.

4.4 Operation

To avoid the risk of personal injury and damage to the column during operation, follow the directions below:

- The pressure specification is valid only if the glass tube is undamaged and free from scratches
- Air should NEVER be used to control the movement of the hydraulic adaptor
- There should be no air or gas in the column

4.5 Maintenance

The O-rings must be intact otherwise chemicals can leak from the column and cause personal injury.

4.6 Hygiene

The INdEX column has not been tested for its biological reactivity according to USP XXIII.

5 Return of equipment to GE Healthcare

For return of equipment's to GE Healthcare for service, maintenance or other reasons please complete enclosed Decontamination Report.

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Decontamination Report

1. Equipment returned for service (or by any other reason), when personnel connected to GE Healthcare should handle the equipment (at any location), should be cleaned.
2. This Form must be used to log the decontamination of the equipment.
3. A completed copy of this form should be Faxed or sent by first class post to the person who will come in contact with the equipment to ensure the he/she have the information before handling the equipment.
4. Failure to complete the form or comply with the procedure will lead to costs for decontamination and delays in handling of the equipment.

5. Company: Address:
Phone: Fax number:

6. Please complete the following sections

6.1 Equipment type 6.2 Serial number

6.3 Details of substances use inside/outside the equipment

6.3.1 Substance/s name:

.....
.....

6.3.2 Precautions to be taken in handling these substances:

.....
.....

6.3.3 Action to be taken in the event of human contact:

.....
.....

6.3.4 Cleaning fluid to be used if residue of chemicals is found during handling:

.....
.....

6.4 Any further relevant information:

.....

7. If substances used in the system are hazardous or toxic, these must be highlighted under 6.3.1.

I hereby confirm that the equipment has been thoroughly cleaned, disregarding if equipment has been in contact with toxic and hazardous substances, and that the equipment has been filled with suitable solution.

I have also attached a description of the cleaning method.

Storage solution in equipment:

Signed: Name:

Position: Date:

www.gehealthcare.com

GE Healthcare Bio-Sciences AB
Björkgatan 30
751 84 Uppsala
Sweden

INdEX is a trademark of GE Healthcare companies. GE, imagination at work and GE monogram are trademarks of General Electric Company.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. GE Healthcare reserves the right, subject to any regulatory and contractual approval, if required, to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your local GE Healthcare representative for the most current information.

© 2006 General Electric Company – All rights reserved.

GE Healthcare Bio-Sciences AB, a General Electric Company.

GE Healthcare Bio-Sciences AB
Björkgatan 30, SE-751 84 Uppsala, Sweden

GE Healthcare Europe GmbH
Munzinger Strasse 5, D-79111 Freiburg, Germany

GE Healthcare UK Ltd
Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK

GE Healthcare Bio-Sciences Corp
800 Centennial Avenue, P.O. Box 1327, Piscataway, NJ 08855-1327, USA

GE Healthcare Bio-Sciences KK
Sanken Bldg. 3-25-1, Hyakunincho, Shinjuku-ku, Tokyo 169-0073, Japan

Asia Pacific Tel: +852 2811 8693 Fax: +852 2811 5251 • Australasia Tel: +61 2 9899 0999 Fax: +61 2 9899 7511 • Austria Tel: 01/57606-1619 Fax: 01/57606-1627 • Belgium Tel: 0800 73 888 Fax: 03 272 1637 • Canada Tel: 800 463 5800 Fax: 800 567 1008 • Central, East, & South East Europe Tel: +43 1 982 3826 Fax: +43 1 985 8327 • Denmark Tel: 45 16 2400 Fax: 45 16 2424 • Finland & Baltics Tel: +358 (0)9 512 39 40 Fax: +358 (0)9 512 39 439 • France Tel: 01 69 35 67 00 Fax: 01 69 41 96 77 • Germany Tel: 0761/4903-490 Fax: 0761/4903-405 • Italy Tel: 02 27322 1 Fax: 02 27302 212 Japan Tel: +81 3 5331 9336 Fax: +81 3 5331 9370 • Latin America Tel: +55 11 3933 7300 Fax: +55 11 3933 7304 • Middle East & Africa Tel: +30 210 9600 687 Fax: +30 210 9600 693 • Netherlands Tel: 0165 580 410 Fax: 0165 580 401 • Norway Tel: 815 65 555 Fax: 815 65 666 • Portugal Tel: 21 417 7035 Fax: 21 417 3184 • Russia & other C.I.S. & N.I.S Tel: +7 (095) 232 0250, 956 1137 Fax: +7 (095) 230 6377 • South East Asia Tel: 60 3 8024 2080 Fax: 60 3 8024 2090 • Spain Tel: 93 594 49 50 Fax: 93 594 49 55 • Sweden Tel: 018 612 1900 Fax: 018 612 1910 • Switzerland Tel: 0848 8028 12 Fax: 0848 8028 13 UK Tel: 0800 616928 Fax: 0800 616927 • USA Tel: 800 526 3593 Fax: 877 295 8102

GE imagination at work

56-3032-64 AD 03/2006